

# Pre-analytical variation in protein biomarker research

## Introduction

Proteins are a critically important class of biomolecules, since they directly control and regulate most of the body's biological functions. They are crucial for health and disease, and the drugs we currently have available largely exert their effects via interactions with proteins. Unlike the genome, the proteome is very dynamic and changes significantly in response to biological signals, environmental conditions and other external stimuli. Given this central role, proteins can serve as invaluable biomarkers to better understand biological mechanisms and to monitor health, therapy effects or the real-time status of any given disease.

Compared to genes, which have been extensively used as biomarkers, proteins are less stable than DNA and therefore require more care in their storage and handling. Taking an extreme example, the DNA of a dinosaur can provide sequence data millions of years after it would have been possible to carry out any reliable analysis of its proteins.



**DNA vs proteins.** You can sequence the DNA in food long after you can't eat it anymore (but storage at the right temperature helps with the latter).

Olink Proteomics has developed high-quality immunoassays for protein biomarker research that can measure many different proteins simultaneously. To ensure the best possible data, however, this assay performance

needs to be matched by good quality, properly handled samples. Proteins can be affected by several factors related to collection and storage, all of which need to be considered when measuring proteins in biological samples and interpreting results. It is important to document as much information as possible about the sample and how it has been collected. Taking this approach, it is perfectly possible to obtain meaningful results even with samples that have not been handled optimally. Proper documentation around sample handling can also provide valuable input when interpreting the data.

This white paper will provide some guidance on what to consider with respect to sample collection and handling when setting up a study. This advice is generally applicable to protein studies using any immunoassay-based approach and is not specific for Olink's technology.

## Serum and plasma

Olink panels can be used for a wide range of sample matrices (see [www.olink.com/sample-types](http://www.olink.com/sample-types) for details), but this document will focus on blood samples. Plasma and serum are both derived from whole blood but a difference in measured levels can be expected for many proteins.

### *Serum*

Serum is the liquid part of blood after it has been allowed to coagulate fully for 30-60 minutes at room temperature. Serum is free of clotting proteins but contains the clotting metabolites that result from the clotting process. The clot is removed by centrifugation.

### *Plasma*

Plasma is the liquid part of blood that has been treated with anti-coagulants, after cells have been removed by centrifugation. Since plasma has been prevented from clotting it is reflective of the blood as it circulates in the body. Plasma collection tubes contain different anti-coagulants such as EDTA, heparin, or sodium citrate, and any of these additives can be used in Olink's analysis, but for consistency, all samples within a study should be treated the same way.



cycles. No significant change among the measured proteins was observed after eight freezing cycles, which indicates that protein levels are stable, and samples with a history of multiple freeze-thaw cycles can still be processed with high quality (1).

Given the challenges of data collection and sample storage within particular studies, there has been little standardization across biobanks. A full data trail on each sample should be provided by the biobank to enable the samples to contribute as part of wider collaborative efforts with other similar samples.

### Transportation

Potential temperature issues arise in the transportation of samples within and between facilities. Care must be taken to ensure that the serum and plasma samples have sufficient dry ice for the expected duration of the transport.

### Future developments

Liquid samples are still the most common type for biomarker research, but filter paper samples with for example dried blood spots, plasma, tear fluids or vaginal fluid are attractive because of the ease and low cost of collection and storage. One study performed by Berglund M. et al. demonstrated the feasibility of measuring proteins from vaginal fluid and plasma dried on filter papers (3). Another study, performed by Björkesten et al., investigated the suitability of dried blood spots, stored at -24°C, for protein measurements. The main findings from the latter were that the act of drying only slightly influenced detection of blood proteins, even after storage for 30 years (4). Dried blood spot biobanks could prove of great medical value by enhancing discovery as well as routine analysis of blood biomarkers.

## How Olink can help

Even with the best intentions and preparations, variations can occur pre-analytically, resulting in individual outliers or overall drift. Olink's technical support can guide you on these matters, and our data science team can help with data analysis. The more information you can provide regarding sample collection, preparation and storage, the better the chances for our data scientists to identify

potential problems to be able to normalize data between samples with different pre-analytical histories.

## Best practice

- *Newly collected plasma samples:* Centrifugate plasma samples as soon as possible, but at least within 1 hour at room temperature or within 8 hours if the samples are kept at 4°C
- *Newly collected serum samples:* Allow serum to fully clot for 30-60 minutes at room temperature prior to centrifugation
- Randomize samples and always analyze the same kind of sample types as cases and controls
- Use dry ice for transportation
- Store samples in a -80°C freezer, as recommended by the global biobanking organization ISBER (5)
- Document how the samples have been collected and handled

## References

1. [Shen et al., Strong impact on plasma protein profiles by precentrifugation delay but not by repeated freeze-thaw cycles, as analyzed using multiplex proximity extension assays. \*Clinical Chemistry and Laboratory Medicine\* \(2017\).](#)
2. [Enroth S., et al., Effects of Long-Term Storage Time and Original Sampling Month on Biobank Plasma Protein Concentrations, \*EBioMedicine\* \(2016\).](#)
3. [Berglund M. et al., Protein Detection Using the Multiplexed Proximity Extension Assay \(PEA\) from Plasma and Vaginal Fluid Applied to the Indicating FTA Elute Micro Card™. \*Journal of Circulation Biomarkers\* \(2016\).](#)
4. [Björkesten et al., Stability of Proteins in Dried Blood Spot Biobanks. \*Molecular & Cellular Proteomics\* 16.7 \(2017\).](#)
5. [International Society for Biological and Environmental Repositories, \*ISBER Best Practices: Recommendations for Repositories Fourth Edition\* \(2018\).](#)

For research use only. Not for use in diagnostic procedures.

This product includes a license for non-commercial use. Commercial users may require additional licenses. Please contact Olink Proteomics AB for details. There are no warranties, expressed or implied, which extend beyond this description. Olink Proteomics AB is not liable for property damage, personal injury, or economic loss caused by this product.

Olink® is a registered trademark of Olink Proteomics AB.

© 2018 Olink Proteomics AB. All third party trademarks are the property of their respective owners.