




# Olink<sup>®</sup> User manual

v 2.0



# Table of contents

1. INTRODUCTION .....	3
About this manual .....	3
Important user information .....	3
Associated documentation .....	3
Technical support .....	4
2. TECHNOLOGY DESCRIPTION .....	5
About PEA .....	5
Quality control .....	5
Data analysis .....	7
3. PRODUCT DESCRIPTION .....	8
Reagents supplied in Olink kits .....	8
4. ADDITIONAL REQUIREMENTS .....	10
NPX manager .....	10
Analysis software .....	10
Required consumables (not supplied) .....	10
Required equipment (not supplied) .....	10
5. LABORATORY WORK SAFETY AND GUIDELINES .....	12
Safety considerations .....	12
Pipetting guidelines .....	12
Vortexing guidelines .....	14
Sterile lab environment .....	15
6. PREPARATION .....	16
Experimental design .....	16
Create programs .....	16
Bench setup .....	17
Time indications and limits .....	17
Sample preparation. ....	18
Plate layout .....	18
7. DILUTION STEP .....	19
Overview .....	19
Sample dilution step for 1:10 panels .....	19
Sample dilution step for 1:100 panels .....	21
Sample dilution step for 1:2025 panels .....	22
8. OPERATION .....	24
Overview .....	24
Incubation step .....	24
Extension and amplification step .....	26
Detection step .....	28
9. OTHER RESOURCES .....	33
Olink FAQs .....	33
Data analysis troubleshooting .....	33
Lab instruction video .....	33
PEA technology video .....	33

# 1. Introduction

Olink high-multiplex immunoassay panels provide an efficient and innovative tool for targeted human protein biomarker discovery, development and validation.

The disease- or biology-focused panels enable rapid, high-throughput analysis, with exceptional data quality and minimal sample consumption. Using just 1 µl of sample, 92 biomarkers are assayed simultaneously on each panel, with results in only 24 hours.

## About this manual

The Olink User Manual provides you with the instructions needed to run an Olink panel.

## Important user information

### *Sample preparation*

To reduce sample-handling time, samples should be distributed in 8-well strips or on a 96-well plate prior to the start of the run.

### *Sample types*

Olink panels have been validated using EDTA plasma and serum samples. A range of additional sample types are compatible with the technology. For example citrate plasma, heparin plasma, tissue and cell lysates, CSF and saliva. Different sample matrices are expected to affect the detection of specific proteins in different ways. In addition, extreme levels of IgG or fluorescent particles can interfere with the Olink assay. For more information on sample types, please see the Data Validation documents corresponding to each panel, or contact Olink support at [support@olink.com](mailto:support@olink.com).

### *Pipettes*

A multichannel pipette and a reverse pipetting technique must be used in the reagent transfer step (see "Pipetting techniques" on page 12). Maintain and calibrate the pipettes regularly.

## Associated documentation

### *Olink documentation*

The *Olink NPX Manager User Guide* is available from the Olink website: [www.olink.com/npx](http://www.olink.com/npx).

### *Other documents*

For information on the Fluidigm IFC Controller HX and Fluidigm BioMark System, read the following User Guides that can be found on [www.fluidigm.com](http://www.fluidigm.com):

- Fluidigm® IFC Controller User Guide - PN 68000112
- Fluidigm® Real-Time PCR Analysis User Guide - PN 68000088
- Fluidigm® Data Collection Software User Guide - PN 68000127


## Technical support

For technical support, please contact Olink Proteomics at [support@olink.com](mailto:support@olink.com).

## 2. Technology description

### About PEA

The technology behind the Olink panels is called Proximity Extension Assay (PEA). The three core steps of the technology are illustrated below.



#### 1. Immuno reaction (Incubation)

Duration: Overnight 16-22 hours

The 92 antibody pairs, labelled with DNA oligonucleotides, bind to their respective protein in the samples.

#### 2. Extension and amplification

Duration: 2 hours

Oligonucleotides that are brought into proximity hybridize, and are extended using a DNA polymerase. This newly created piece of DNA barcode is amplified by PCR.

#### 3. Detection

Duration: 4.5 hours


The amount of each DNA barcode is quantified by microfluidic qPCR.

### Quality control

Olink has developed a built-in QC system using internal controls, for its multiplex biomarker panels. This system enables full control over the technical performance of assays and samples.

### Internal controls

The QC system consists of four internal controls that are spiked into every sample and designed to monitor the three main steps of the Olink protocol: Immuno reaction (incubation step), extension and amplification/detection.




**Incubation controls:** Incubation Control 1 and 2 are two different non-human antigens measured with PEA. These controls monitor potential technical variation in all three steps of the reaction.

**Extension control:** The Extension Control is composed of an antibody coupled to a unique pair of DNA-tags. These DNA-tags are always in proximity, so that this control is expected to give a constant signal independently of the immuno reaction. This control monitors variation in the extension and amplification/detection step and is used to adjust the signal from each sample with respect to extension and amplification.

**Detection control:** The Detection Control is a complete double stranded DNA amplicon which does not require any proximity binding or extension step to generate a signal. This control monitors the amplification/detection step.

### External controls


There are six required and two recommended external controls that are added to separate wells on each sample plate. On the illustrated plate there are 88 samples, 2 pooled plasma samples called External sample controls, 3 Negative controls and 3 Inter-plate controls.



**Inter-plate Control:** Inter-plate Control (IPC) is included in triplicate on each plate and these are run as normal samples. The IPC is a pool of 92 antibodies, each with one pair of unique DNA-tags positioned in fixed proximity and can be seen as a synthetic sample, expected to give a high signal for all assays. The median of the IPC triplicates is used to normalize each assay, to compensate for potential variation between runs and plates.

**Negative Control:** Negative Control is also included in triplicate on each plate and consists of buffer run as a normal sample. These are used to monitor any background noise generated when DNA-tags come in close proximity without prior binding to the appropriate protein. The negative controls set the background levels for each protein assay and are used to calculate the limit of detection (LOD).

**Sample Control:** On each plate, it is recommended to run a pooled plasma sample in duplicate. These are used to assess potential variation between runs and plates, for example to calculate inter-assay and intra-assay CV, as well as troubleshooting. Validation data for all panels is available on the Olink website. An example of inter- and intra-assay CV is illustrated here:



## Data analysis

### NPX values

Olink reports protein expression levels using an arbitrary unit called Normalized Protein eXpression (NPX). NPX is a relative quantification unit on log<sub>2</sub> scale that allows users to identify changes for individual protein levels across their sample set, and then use this data to establish protein signatures. NPX is derived from the Ct values obtained from the qPCR using the following equations:

#### Extension Control:

$$Ct_{Analyte} - Ct_{Extension\ Control} = dCt_{Analyte}$$

#### Inter-plate Control:

$$dCt_{Analyte} - dCt_{Inter\ plate\ Control} = ddCt_{Analyte}$$

#### Adjustment against a correction factor:

$$Correction\ factor - ddCt_{Analyte} = NPX_{Analyte}$$

The data processing (normalization) of the Ct-values is performed to minimize both intra- and inter-assay variation and make the data more intuitive and easy to interpret. Quality control and normalization is achieved using the Olink NPX Manager software (see "NPX manager" on page 10).

The NPX unit is unique to each protein assay, meaning that even if two different proteins have the same NPX values, their absolute concentrations may differ. Due to the relative nature of the unit, NPX should not be compared between runs without proper inter-plate normalization

See the Olink white paper *Data normalization and standardization* for more information (download from [www.olink.com/whitepapers](http://www.olink.com/whitepapers)).

# 3. Product description

## Reagents supplied in Olink kits

Each Olink kit contains reagents for 96 wells, sufficient for 90 samples and 6 controls. The reagents are supplied in two individual boxes. Storage temperature and expiry date for the components are stated on the label on each box.



### Small box

The small box includes the Olink Probe kit and should be stored at +4°C. The Probe kit includes:

Part	Description
Incubation Solution	Contains components needed for the incubation reaction
A-probes	Contains 96 antibody probes labeled with A oligos
B-probes	Contains 96 antibody probes labeled with B oligos

### Large box

The large box contains the Olink Detection and Control kit. It should be stored at -20°C and includes:

Part	Description
PEA Solution	Contains components needed for the extension reaction
PEA Enzyme	For extension of A and B probes bound to their target
PCR Polymerase	For pre-amplification of the extension product created by the PEA Enzyme, also used in the Detection step
Detection Solution	Contains components needed for the detection reaction
Detection Enzyme	For qPCR amplification
Primer Plate	96-well plate with ready-to-use primers for amplification of extension product
Inter-plate Control	For normalization of each assay, to compensate for potential variation between runs and plates
Negative Control	For determination of background levels
Incubation Stabilizer	For stabilization of the incubation reaction
Sample Diluent	Only included for panels that require pre-dilution of samples

## 4. Additional requirements

### NPX manager

Olink NPX Manager software is an easy to use data import and pre-processing tool developed by Olink Proteomics. The software lets you import data, validate data quality and normalize Olink data for subsequent statistical analysis. See the *Olink NPX Manager User Guide* for more information.

### Analysis software

The Fluidigm® Biomark™ Data Collection software, RT-PCR is required for the analysis of Olink data in the NPX Manager software.

#### *Protein list*

The list of proteins can directly be imported into the Fluidigm analysis software as a .plt file. The relevant .plt file can be found on the lot configuration data sheet provided in every kit. Download the Detector Setup.plt file corresponding to your panel at [www.olink.com/data](http://www.olink.com/data).

### Required consumables (not supplied)

**Tip:** Contact Olink support at [support@olink.com](mailto:support@olink.com) for specific recommendations.

- Pipette filter tips
- Microcentrifuge tubes (1-1.5 mL)
- Centrifuge tube (>11 mL)
- 8-well strips with lids
- 96-well PCR plate (à 0.2 mL)
- Multichannel pipette reservoir
- Adhesive plastic film (heat-resistant)
- High purity water (sterile filtered, MilliQ® or similar)
- 96.96 Dynamic Array™ Integrated Fluidic Circuit (IFC), (Fluidigm Corporation, catalogue number BMK-M96.96)
- External Sample Controls (Pooled plasma samples)

### Required equipment (not supplied)

- Pipettes (covering the range from 1 µL to 1000 µL) and compatible pipette tips
- Multichannel pipettes (recommended range 1-10 µL and 50-100 µL and/or 50-200 µL) and compatible multichannel pipette tips
- Vortex

- Centrifuge for plates
- Microcentrifuge for tubes
- Freezing block (-20°C) for enzyme handling
- Thermal cycler with:
  - Heated lid
  - Temperature range from +50°C to +95°C
  - Validated for 0.1 mL volumes (important)
  - 96-well format (recommended)
- Refrigerator or cold room (+2°C to +8°C)
- Freezer (-20°C)
- Fluidigm BioMark™ or BioMark™ HD System
- Fluidigm IFC controller HX

# 5. Laboratory work safety and guidelines

## Safety considerations


### Safety

Follow general laboratory safety procedures:

- Use gloves, safety goggles and protective clothing when performing the experiments.
- Handle and dispose of hazardous sample material according to local regulations.

### PCR technology

PCR technology is sensitive to contaminations. Perform the Detection step in a PCR room, separate from the room where the Incubation and Extension steps are performed. Maintain and calibrate the PCR and BioMark™ instruments regularly.



## Pipetting guidelines

### Tips and tricks for pipetting

- Calibrate pipettes regularly (at least a 6 month interval)
- Pipette near the liquid surface
- Let the reagents and liquids reach room temperature before use to maximize accuracy.
- Pre-rinsing the tip 1 to 3 times with the liquid to be pipetted improves accuracy when using reverse pipetting
- Do not turn the pipette on its side when there is liquid in the tip as liquid might contaminate the interior of the pipette
- Keep pipettes vertical while pipetting and pipette to the bottom of the wells.

### Pipetting techniques

Both forward and reverse pipetting is used in the Olink protocol.

Forward pipetting is the most commonly used pipetting technique and is performed as described in the following instruction:

1. Press the operating button to the first stop.

- Dip the tip into the solution to a depth in accordance with the set volume, and slowly release the operating button. Remove the tip from the liquid.

**NOTE:** Ensure that all tips contain the exact same volume if a multichannel pipette is used.

- Dispense the liquid into the receiving vessel by gently pressing the operating button to the first stop and then press the operating button to the second stop. This action will empty the tip. Remove the tip from the vessel.
- Release the operating button to the ready position.

Ready position	1	2	3	4
First stop	↓	↑	↓	↑
Second step			↓	↑

Reverse pipetting improves precision with smaller volumes and viscous solutions and is performed as follows:

- Press the operating button past the first stop.
- Dip the tip into the solution to a depth in accordance with the volume set, and slowly release the operating button. This action will fill the tip with a volume that is larger than the set volume.

**NOTE:** Ensure that all tips contain the exact same volume if a multichannel pipette is used.

- Remove the tip from the liquid and dispense the liquid into the receiving vessel by pressing the operating button gently and steadily to the first stop. This volume is equal to the set volume.
- Hold the button in this position. Some liquid will remain in the tip- this should not be dispensed.
- Continue pipetting by repeating steps 3 and 4.

Ready position	1	2	3	4	5	...	X	End
First stop	↓	↑	↓	↑	↓			↑
Second step	↓	↑					↓	↑

## Vortexing guidelines

Correct vortexing is critical when running Olink panels.


### *Vortexing tips and tricks*

- Cover the rubber platform of the vortex with adhesive plastic film to make it easier to slide the 96-well plate during vortexing
- Vortex for 20-30 seconds at full speed
- Visually inspect the wells during vortexing to ensure complete mixing. The liquid should swirl in the wells.


### *Vortexing instruction*

Correct vortexing is essential for generating reproducible results. To vortex thoroughly, follow these steps:


1. Move the plate back and forth over the vortex in horizontal direction. Make sure that the wells at the outside edges of the plate (columns 1 and 12) are also vortexed.



2. Tilt the plate away from you and move it back and forth over the vortex.



3. Tilt the plate towards you and move it back and forth over the vortex.



4. Turn the plate 180° and repeat steps 1-3.

## Sterile lab environment

The high sensitivity of the Olink assays requires a clean laboratory environment. Particles from the surroundings, such as dust, hair, saliva and skin flakes are common sources of contamination. The following recommendations are intended to simplify the workflow in the lab and improve data quality.

- Keep all consumables (tubes, pipette tips, PCR plates etc.) in closed bags or boxes, preferably in a closed storage unit.
- Always wear gloves, including when bringing reagents in and out from fridge or freezer.
- Change gloves when needed and always when going from pre-PCR to PCR operations.
- Always wear a long-sleeved lab coat.
- Use a separate lab coat for working with the PCR steps.
- Use separate rooms for pre-PCR and PCR operations. If this is not possible, keep separate benches.
- Use separate consumables and equipment for pre-PCR and PCR operations.

# 6. Preparation

## Experimental design

Decide how many samples, replicates and controls that are needed to get the data you want from the study. When running more than one plate it is important that the samples are appropriately randomized across all plates necessary steps for normalizing and combining data are taken. It may be wise to consult a statistician or Olink Support ([support@olink.com](mailto:support@olink.com)) prior to running the study. See the Olink white paper *Strategies for design of protein biomarker studies* for more information (download from [www.olink.com/whitepapers](http://www.olink.com/whitepapers)).

## Create programs

### *Olink extension and pre-amplification program*

Create a PEA program on the thermal cycler with the following settings:

Step	Temperature	Duration
Extension	50°C	20 min
Hot start	95°C	5 min
PCR cycle	95°C	30 s
	54°C	1 min
	60°C	1 min
Maintain the reaction at	10°C	∞, hold

} ×17

**NOTE:** Enable the heated lid function.

### *Olink Protein Expression 96x96 v2 program*

Create a Fluidigm BioMark System program named "Olink Protein Expression 96x96 v2" with the following steps:

Step	Temperature	Duration
Thermal mix	50°C	120 s
	70°C	1800 s
	25°C	600 s
Hot start	95°C	300 s
PCR cycle	95°C	15 s
	60°C	60 s

} ×35

Verify that the program has the following settings:

Variable	Setting
Application	Gene Expression
Passive Reference	5-Carboxy-x-Rhodamine (abbreviation ROX in Fluidigm software)
Assay	Single probe
Probes	FAM-MGB

## Bench setup

Before you start working it is important to know that you have everything you need and that all reagents are ready for pipetting. Do the following before you start the experiment:

1. Clean the lab bench, hood, racks and pipettes with 70% ethanol.
2. Bring out all reagents, consumables and samples needed for the specific lab step. Potential contamination from the surrounding environment is minimized by reducing physical movements in the lab.
3. Organize equipment, consumables and samples in the work station, in a way that enables clean work.
4. Label pipette boxes with column number to more easily monitor where you are on the plate.
5. Note that all consumables and reagents are single use only.

## Time indications and limits

Step	Duration	Comment
Dilution (for diluted panels)	15-30 min	The diluted panels Metabolism, CVD III, Development and Cardiometabolic need a dilution step.
Incubation setup	30-45 min	
Incubation reaction	16-22 hours (overnight)	Keep incubation times consistent when running multiple chips for the same project (variation <2 hours).
Extension	30 min	From preparing mix to start of PCR machine <b>NOTE:</b> A maximum of 5 minutes can pass between adding extension mix to the incubation plate until the start of the PCR. This is important for sensitivity.
Priming	20 min	Prepare sample plate during priming of the chip
Load	2 hours	Start adding primers and samples to primed chip within 10 minutes of completed priming Start the load script in the IFC controller within 40 minutes The time in the IFC controller is 90 minutes.
qPCR	2.5 hours	Move the loaded chip from the IFC controller to the Biomark within 30 minutes of loading completion <b>NOTE:</b> We do not recommend setting up more than 3 plates at a time.


## Sample preparation.

Sample preparation recommendations are listed below:

- Use a 96-well PCR plate format, preferably with full skirt. All plates need to be able to withstand -80°C, be dry-ice resistant and easily re-sealable.
- Ensure that the samples are randomized, or in the order ready to be run and compatible with the plate layout.
- Ensure that each well is separately sealed using an adhesive film or individual seals.
- Clearly mark sample plates or tubes with a simple alphanumeric code that you can later identify (“A, B, C”, “1, 2, 3”, or “A1, A2, A3”) using temperature-resistant labels or marker pen.
- Use unique sample identification names or numbers.

## Plate layout

Below is an illustration of the plate layout. The twelfth column consists of controls. The negative controls and inter-plate controls are always supplied by Olink. The pooled sample controls are provided by Olink if the samples are run at an Olink Analysis Service lab.



● Negative Control ● Inter-plate Control ● Sample Control ● Samples

# 7. Dilution step

## Overview

For selected panels, a sample dilution step is required prior to running the assay protocol. For Olink panels that require pre-dilution of samples, this information is indicated on the kit boxes as well as on the lot configuration sheet. During this step, the samples are diluted so that the target proteins are in an optimal concentration range for the assay.

**NOTE:** Dilutions are optimized for serum and plasma only. Other sample matrices may require different dilutions. Contact [support@olink.com](mailto:support@olink.com) for more information.

**NOTE:** Always double check that you have the correct volume in ALL pipette tips!

Dilution	Panels
1:10	Metabolism
1:100	CVD III and Development
1:2025	Cardiometabolic

## Sample dilution step for 1:10 panels

### Prepare bench

For the Dilution step for 1:10 panels, you will need:

- Prepared 96-well plate with samples
- 1 96-well PCR plate
- 1 Sample Diluent
- 1 multichannel pipette reservoir
- 1 multichannel pipette (10  $\mu$ L)
- 98 multichannel pipette tips (10  $\mu$ L)
- 2 adhesive films

### Before you start

- Decide how many samples you will include in the experiment.
- Use the 96-well plate layout and select a location for each sample. Make sure the samples are distributed randomly.

**NOTE:** The negative control and inter-plate control samples should **not** be diluted.

- Sample dilutions should be made in a 96-well plate (0.2 mL per well) using a multichannel pipette.

**NOTE:** *Pipette near the surface of the samples to prevent liquid from sticking to the outside of the pipette tip.*

- Prepare the dilutions freshly, as close to the start time of the assay as possible, no more than 1 hour before the start of incubation. Store the dilutions in the fridge.

### *Dilution step instruction*

1. Thaw the Sample Diluent, vortex and empty the bottle into a multichannel pipette reservoir (minimum volume 15 mL).
2. Mark a 96-well plate as "Dilution Plate".
3. Pipette up and down a few times to pre-condition the pipette tips.
4. Transfer 9  $\mu$ L of the Sample Diluent to each well of columns 1-11 and positions A-B in column 12 on the 96-well plate, using **reverse pipetting**. Pipette the Sample Diluent carefully to avoid foaming.

**NOTE:** *Do not change pipette tips and use the same multichannel pipette throughout the entire dilution series.*

5. Vortex the samples and spin down the liquid at 400 x g, for 1 minute at room temperature. Carefully transfer 1  $\mu$ L of your samples and pooled sample controls according to your plate layout to the Dilution Plate using **forward pipetting**.

**NOTE:** *Use the same multichannel pipette throughout the entire plate, also for samples A-B of column 12. Change tips between each pipetting step.*

6. Seal both the original sample plate and the Dilution Plate with adhesive plastic film.
7. Vortex the Dilution Plate thoroughly. Move the plate on the vortex and tilt the plate to make sure that the liquid is moving in the wells. See "Vortexing guidelines" on page 14 for more information.
8. Spin down the content at 400 x g for 1 minute at room temperature.
9. Double check that all wells in the plate contain the same volume. Note any deviations. This information is needed to interpret the data in the statistical analysis step.

## Sample dilution step for 1:100 panels

### *Prepare bench*

For the Dilution step for 1:100 panels, you will need:

- Prepared 96-well plate with samples
- 1 96-well PCR plate
- 1 Sample Diluent
- 1 multichannel pipette reservoir

- 1 multichannel pipette (200  $\mu$ L)
- 1 multichannel pipette (10  $\mu$ L)
- 8 multichannel pipette tips (200  $\mu$ L)
- 90 multichannel pipette tips (10  $\mu$ L)
- 2 adhesive films

### *Before you start*

- Decide how many samples you will include in the experiment.
- Use the 96-well plate layout and select a location for each sample. Make sure the samples are distributed randomly.

**NOTE:** *The negative control and inter-plate control samples should **not** be diluted.*

- Sample dilutions should be made in a 96-well plate (0.2 mL per well) using a multichannel pipette.

**NOTE:** *Pipette near the surface of the samples to prevent liquid from sticking to the outside of the pipette tip.*

- Prepare the dilutions freshly, as close to the start time of the assay as possible, no more than 1 hour before the start of incubation. Store the dilutions in the fridge.

### *Dilution step instruction*

1. Thaw the Sample Diluent, vortex and empty the bottle into a multichannel pipette reservoir (minimum volume of 15 mL).
2. Mark a 96-well plate as "Dilution Plate".
3. Pipette up and down a few times to pre-condition the pipette tips.
4. Transfer 99  $\mu$ L of the Sample Diluent to each well of columns 1-11 and positions A-B in column 12 on the 96-well plate, by using **reverse pipetting**. Pipette the Sample Diluent carefully to avoid foaming.

**NOTE:** *Do not change pipette tips and use the same multichannel pipette throughout the entire dilution series.*

5. Vortex the samples and spin down the liquid at 400 x g, for 1 minute at room temperature. Carefully transfer 1  $\mu$ L of your samples and pooled sample controls according to your plate layout using **forward pipetting**.

**NOTE:** *Use the same multichannel pipette throughout the entire plate, also for samples A-B of column 12. Change tips between each pipetting step.*

6. Seal both the original sample plate and the Dilution Plate with adhesive plastic film.
7. Vortex the Dilution Plate plate thoroughly. Move the plate on the vortex and tilt the plate to see that the liquid is moving in the wells. See "Vortexing guidelines" on page 14 for more information.
8. Spin down the content at 400 x g for 1 minute at room temperature.

9. Double check that all wells in the plate contain the same volume. Note any deviations. This information is needed to interpret the data in the statistical analysis step.

## Sample dilution step for 1:2025 panels

### *Prepare bench*

For the Dilution step for 1:2025 panels, you will need:

- Prepared 96-well plate with samples
- 2 96-well PCR plates
- 1 Sample Diluent
- 1 multichannel pipette reservoir
- 1 multichannel pipette (200 µL)
- 1 multichannel pipette (10 µL)
- 16 multichannel pipette tips (200 µL)
- 180 multichannel pipette tips (10 µL)
- 3 adhesive films

### *Before you start*

- Decide how many samples you will include in the experiment.
- Use the 96-well plate layout and select a location for each sample. Make sure the samples are distributed randomly.

**NOTE:** *The negative control and inter-plate control samples should **not** be diluted.*

- A 2-step sample dilution should be performed in order to reach the desired dilution factor of 1:2025, and made in 2×96-well plates (0.2 mL per well) using a multichannel pipette.

**NOTE:** *Pipette near the surface of the samples to prevent liquid from sticking to the outside of the pipette tips.*

- Prepare the dilutions freshly, as close to the start time of the assay as possible, no more than 1 hour before the start of incubation. Store the dilutions in the fridge.

### *Dilution step instruction*

1. Thaw the Sample Diluent, vortex and empty the bottle into a multichannel pipette reservoir (minimum volume of 15 mL).
2. Mark two 96-well plates as "Dilution Plate 1" and "Dilution Plate 2".
3. Pipette up and down a few times to pre-condition the pipette tips.
4. Transfer 44 µL of the Sample Diluent into each well of columns 1–11 and positions A–B in column 12 on Dilution Plate 1 and Dilution Plate 2 using **reverse pipetting**.

**NOTE:** Do not change pipette tips and use the same multichannel pipette throughout the entire dilution series.

**NOTE:** Pipette the Sample Diluent carefully to avoid foaming.

5. Vortex the samples and spin down the liquid at 400 x g, for 1 minute at room temperature. Carefully transfer 1  $\mu$ L of your samples and pooled sample controls into Dilution Plate 1 according to your plate layout using **forward pipetting**. See "Plate layout" on page 18 for more information.

**NOTE:** Use the same multichannel pipette throughout the entire plate, also for samples A-B of column 12. Change tips between each pipetting step.

6. Seal both the original sample plate and Dilution Plate 1 with an adhesive plastic film.
7. Vortex Dilution Plate 1 thoroughly and ensure that all wells are mixed. Spin down the contents at 400  $\times$  g for 1 minute at room temperature.
8. Transfer 1  $\mu$ L of your diluted samples and pooled sample controls from Dilution Plate 1 to the previously prepared Dilution Plate 2 using **forward pipetting**.


**NOTE:** Use the same multichannel pipette throughout the entire plate, also for samples A-B of column 12. Change tips between each pipetting step.

9. Seal Dilution Plate 2 with adhesive plastic film.
10. Vortex Dilution Plate 2 thoroughly and ensure that all wells are mixed. Spin down the contents at 400  $\times$  g for 1 minute at room temperature.
11. Double check that all wells in the plate contain the same volume. Note any deviations. This information is needed to interpret the data in the statistical analysis step.

# 8. Operation

## Overview

The Olink protocol consists of three core steps. These are the Incubation step, the Extension and amplification step and the Detection step. Panels that target high abundant proteins require an additional Sample dilution step. See "Dilution step" on page 19.



### 1. Incubation

Duration: Overnight 16-22 hours

Location: Pre-PCR room

Needed for this step:

- A-probes
- B-probes
- Incubation Solution
- Incubation Stabilizer
- Negative Controls
- Inter-plate Controls
- Pooled sample controls

### 2. Extension and amplification

Duration: 2 hours

Location: Pre-PCR and PCR room

Needed for this step:

- PEA Solution
- PEA Enzyme
- PCR Polymerase

### 3. Detection

Duration: 4,5 hours

Location: PCR room

Needed for this step:

- Detection Solution
- Detection Enzyme
- PCR Polymerase
- Primer plate
- Biomark chip

## Incubation step

The Incubation step is where the antibody-pairs with attached DNA tags are added to the samples, and allowed to bind to their target proteins during an overnight incubation.

### Prepare bench

Prepare the lab space with pipettes and consumables needed for this step and let the reagents and controls reach room temperature. For the Incubation step, you will need:

- Incubation Solution
- Incubation Stabilizer
- A- and B-probes
- Negative and Inter-plate Controls.
- Pooled sample controls
- 2 8-well strips
- Prepared 96-well plate with samples

- 1 96-well plate
- 1 microcentrifuge tube (1-1.5 mL)
- 3 pipette tips (10  $\mu$ L)
- 4 pipette tips (100/200  $\mu$ L)
- 1 pipette tip (1000  $\mu$ L)
- 104 multichannel pipette tips (10  $\mu$ L)
- 2 adhesive films


**Vortex and spin all reagents before use.**

### *Incubation step instruction*

**Tip:** Practice reverse pipetting before you start. See "Pipetting techniques" on page 12.

**NOTE:** Always double check that you have the correct volume in ALL pipette tips!

1. Thaw the samples.
2. Vortex the Negative Control, Inter-plate Control and Sample Control and spin briefly. Add 5  $\mu$ L of the controls to an 8-well strip according to the following structure:



**NOTE:** If running a diluted panel, the two sample controls will be placed on the sample plate and not in the control strip.

3. Prepare the Incubation Mix in a microcentrifuge tube according to the table below. Vortex and spin each reagent before adding it to the mix.

**NOTE:** Pipette the Incubation Solution carefully to avoid foaming.

Incubation Mix	Per 96-well plate ( $\mu$ L)
Incubation Solution	280.0
Incubation Stabilizer	40.0
A-probes	40.0
B-probes	40.0
<b>Total</b>	<b>400.0</b>

4. Vortex and spin down the Incubation Mix. Transfer 47  $\mu$ L of the Incubation Mix to each well of a new 8-well strip.
5. Pre-condition the pipette tips and transfer 3  $\mu$ L of Incubation Mix to the bottom of the

wells of a new 96-well plate by **reverse pipetting** and name it Incubation Plate. Use the same pipette tips for the entire plate.

**NOTE:** *Pipette near the surface of the Incubation Mix to prevent liquid from sticking to the outside of the pipette tips.*

6. Vortex the samples and spin down the liquid at 400 x g, for 1 minute at room temperature. Transfer 1  $\mu\text{L}$  of each sample, using a multichannel pipette, to the bottom of the wells of the Incubation plate according to your sample plate layout.

**NOTE:** *For diluted panels, transfer from the final Dilution Plate.*

Use **forward pipetting** and change pipette tips between every column.

7. Use a multichannel pipette to transfer 1  $\mu\text{L}$  of Negative Control and Inter-plate Control from the prepared 8-well strip with controls, to the last 6 wells of column 12 of the Incubation Plate, according to the plate layout. See "Plate layout" on page 18 for more information. Use **forward pipetting**.
8. Seal the Incubation Plate thoroughly with an adhesive plastic film, spin at 400 x g for 1 minute at room temperature.

**NOTE:** *Do not vortex the plate.*

9. Double check that all wells in the plate contain the same volume. Note any deviations. This information is needed to interpret the data in the statistical analysis step.
10. Incubate the Incubation Plate overnight at +4°C for 16-22 hours.

**NOTE:** *It is important that all wells are properly sealed to avoid evaporation of the samples.*

## Extension and amplification step

In the morning of the following day, the extension and amplification steps take place. Unique DNA reporter sequences for each target protein are generated and pre-amplified using regular PCR.

### Prepare bench

Prepare the lab space with pipettes and consumables needed for the Extension and amplification step. You will need:

- High purity water
  - PEA Solution
  - PEA Enzyme
  - PCR Polymerase
- NOTE:** *Keep the PEA Enzyme and the PCR Polymerase in a freezing block or on ice.*
- 1 centrifuge tube (>11 mL)
  - 1 multichannel pipette reservoir

- 2 pipette tips (1000 µL)
- 2 pipette tips (100/200 µL)
- 8 multichannel pipette tips (200 µL)
- 1 adhesive film

### Extension and amplification step instruction

**NOTE:** Always double check that you have the correct volume in ALL pipette tips!

1. Allow the PEA Solution to reach room temperature. Vortex and spin down briefly before use.
2. Pre-heat the PCR machine to 50°C and pause the program.
3. Spin down the Incubation Plate at 400 x g for 1 minute at room temperature.
4. Prepare the Extension Mix in a 15 mL tube according to the following table:

Extension Mix	Per 96-well plate (µL)
High Purity Water	9385
PEA Solution	1100
PEA Enzyme	55
PCR Polymerase	22
<b>Total</b>	<b>10 562</b>

**NOTE:** Vortex and spin down all the Extension Mix components tubes.

**Save the remainder of the PCR Polymerase for the final detection step.**

5. Vortex the Extension Mix and pour it into a multichannel pipette reservoir.
6. Carefully remove the adhesive film from the Incubation Plate.
7. Start a 5 minute timer and transfer 96 µL of Extension Mix to the upper parts of each of the well walls of the Incubation Plate using **reverse pipetting**. Use the same pipette tips throughout the plate.

**NOTE:** Do not let the tips come in contact with the content of a well.

8. Seal the plate with a new adhesive plastic film.

**NOTE:** It is important that all wells are properly sealed to avoid evaporation of the samples.

9. Vortex the plate thoroughly to ensure that all wells are mixed before spinning it down. See "Vortexing guidelines" on page 14 for more information.
10. Double check that all wells in the plate contain the same volume. Note any deviations. This information is needed to interpret the data in the statistical analysis step.
11. Take the Incubation Plate to the PCR room.



Perform steps 7-11 within 5 minutes.

12. Immediately place the Incubation Plate in the PCR instrument and resume the PEA program. The PEA program takes approximately 1 hour and 30 minutes. See "Olink extension and pre-amplification program" on page 16 for more information.

**NOTE:** *If the thermal cycler requires a silicon cover for plates covered with plastic film, use one to avoid evaporation.*

**Tip:** *You can start the preparations for the detection step during the last 10 minutes of the PEA program.*

13. When the program is finished, you can continue with the Detection step or store the extension products in the Incubation Plate for up to one week at +4°C.

## Detection step

The final Detection step quantifies the DNA reporters for each biomarker using high throughput real-time qPCR on the Fluidigm Biomark system.

### *Prepare bench*


For the priming of the chip, you need:

- 96.96 Dynamic Array IFC
- Two syringes with control line fluid

For the Detection step, you need:

- Detection Solution
- High purity water
- Detection Enzyme
- PCR Polymerase
- Primer Plate
- Extension products
- 1 8-well strip
- 1 96-well plate
- 1 microcentrifuge tube (1-1.5 mL)
- 2 pipette tips (10 µL)
- 1 pipette tip (100/200 µL)
- 2 pipette tips (1000 µL)
- 3 boxes + 8 multichannel pipette tips (10 µL)
- 2 adhesive films

## Prime chip



1. Unpack the syringes and actuate both check valves of the chip with gentle pressure using a syringe. Ensure that the poppets of the check valves can move freely up and down.
2. Tilt the chip and insert the syringe tip into the check valve opening. The tip should be fully inserted. Confirm that the valve is fully open with the O-ring seal pushed down and moved to the side.
3. Slowly inject one control line fluid syringe into each accumulator on the chip.
4. Place the chip with barcoding facing you in the IFC Controller. First select **Load Chip** and then **Prime** on the IFC Controller. This program takes approximately 20 minutes.

**NOTE:** During the priming of the chip, you should prepare the Detection Mix.

## Detection step instruction

**NOTE:** Always double check that you have the correct volume in ALL pipette tips!

1. Remove the Extension products from the PCR instrument. Vortex and spin down the liquid.
2. Double check that all wells in the plate contain the same volume. Note any deviations. This information is needed to interpret the data in the statistical analysis step.
3. Thaw the Primer Plate and Detection Solution. Vortex and spin down the Detection Solution and Primer Plate. Keep the Detection Enzyme and PCR Polymerase in a freezing block or on ice.


4. Prepare the Detection Mix in a microcentrifuge tube. Vortex and spin down the enzymes before adding them to the mix.

Detection Mix	Per 96-well plate ( $\mu\text{L}$ )
Detection Solution	550.0
High purity water	230.0
Detection Enzyme	7.8
PCR Polymerase	3.1
<b>Total</b>	<b>790.9</b>

5. Vortex the Detection Mix and spin briefly. Transfer 95  $\mu\text{L}$  of the mix to each well of an 8-well strip.
6. Use a multichannel pipette to transfer 7.2  $\mu\text{L}$  of the Detection Mix to each well of a new 96-well plate by **reverse pipetting**. Use the same pipette tips throughout the plate. Name this plate Sample Plate.
7. Carefully remove the adhesive film from the Incubation Plate.
8. Transfer 2.8  $\mu\text{L}$  from the extension products in the Incubation Plate to the Sample Plate using a multichannel pipette and **forward pipetting**. Change tips between each column.
9. Seal both the Sample Plate and the plate with extension products with adhesive plastic film.  
*Tip: The plate with extension products can be saved for 1 week at +4°C.*
10. Vortex and spin the Sample Plate at 400 x g, for 1 minute at room temperature together with the Primer Plate.
11. Double check that all wells in the Sample Plate contain the same volume. Note any deviations. This information is needed to interpret the data in the statistical analysis step.
12. Remove the primed chip from the IFC controller.


**NOTE:** *The chip should be oriented so that the cut corner of the chip is placed on the upper left side.*

In the following steps, the primers will be loaded to the left and samples to the right on the 96.96 Dynamic Array IFC. Here follows an overview illustration of the loading:




13. Carefully remove the adhesive film from the Primer Plate to avoid contamination between wells.

14. Transfer 5  $\mu\text{L}$  using **reverse pipetting** from each well in the Primer Plate to the inlets on the left side of the chip. Change pipette tips after each primer.




15. Transfer 5  $\mu\text{L}$  of each sample using **reverse pipetting** from each well in position 1 A-H (green) to the inlets in the first column on the right side of the chip (green). Change pipette tips after each sample. When using an eight-channel pipette every other inlet will be filled according to the image.




16. Transfer 5  $\mu\text{L}$  from each well in position 2 A-H (blue) to the second column of inlets (blue) according to the image. Continue with columns 3-6. Use **reverse pipetting** and

change tips between each column.



17. Transfer 5  $\mu$ L from each well in position 7 A-H (red) to the inlets in the first column on the right side of the chip (red), start on the second row according to image. Continue with columns 8-12. Use **reverse pipetting** and change tips between each column.



18. When loading is finished, inspect the wells and remove any bubbles using a syringe needle. Change needle between wells to avoid contamination.
19. Place the chip with its barcoding facing you in the IFC Controller. Select **Load Chip** and **Load Mix** followed by **Run Script** to load the assay and sample mixes into the central part of the chip. The loading program takes approximately 90 minutes. See the manufacturer's instructions for more details.
20. When the chip loading is completed, eject the chip from the IFC controller.
21. Use a piece of adhesive tape to remove dust from the top of the middle section of the chip. Let the sticky part lightly touch the surface of the chip.
22. Remove the protective film from underneath the chip.
23. Load the chip in the Fluidigm BioMark with the barcode facing outwards and start the Olink Protein Expression protocol. See "Olink Protein Expression 96x96 v2 program" on page 16. The Biomark program takes approximately 2 hours and 10 minutes.

## 9. Other resources

### Olink FAQs

The answers to the most common questions asked by our customers can be found on the Olink website: [www.olink.com/faq](http://www.olink.com/faq).

### Data analysis troubleshooting

Regarding issues detected during the quality control of the data or statistical analysis, see the troubleshooting chapter in the *Olink NPX Manager User Guide* for solutions. The *Olink NPX Manager User Guide* is available from the Olink website: [www.olink.com/npx](http://www.olink.com/npx).

### Lab instruction video

Learn from our experienced Analysis Service team how to run Olink panels by watching the [Olink Lab Instruction video](#) on the Olink youtube channel.

### PEA technology video

For an animated description of how our innovative dual recognition, DNA-coupled methodology provides exceptional readout specificity, watch the [PEA overview](#) on the Olink youtube channel.

For research use only. Not for use in diagnostic procedures.  
This product includes a license for non-commercial use. Commercial users may require additional licenses. Please contact Olink Proteomics AB for details.  
There are no warranties, expressed or implied, which extend beyond this description. Olink Proteomics AB is not liable for property damage, personal injury, or economic loss caused by this product.  
Olink® is a registered trademark of Olink Proteomics AB.  
Biomark™ and Fluidigm® are trademarks of Fluidigm Corporation.  
© 2019 Olink Proteomics AB. All third party trademarks are the property of their respective owners.



Olink Proteomics  
Dag Hammarskjölds v. 52B  
SE-752 37 Uppsala, Sweden  
[www.olink.com](http://www.olink.com)