

Olink®

COVID-19 protocol certification

Olink Proteomics AB is committed to partnering with the scientific community to learn more about COVID-19 in hopes of finding a therapeutic cure. To that end, Olink's Analysis Service Uppsala and Boston office will accept COVID-19 patient samples for testing provided the requirements set forth below are STRICTLY followed, as certified by an authorized representative of your institution. ALL samples sent to Olink for analysis must undergo the same inactivation protocol, REGARDLESS of COVID-19 status (positive or negative):

- A. Accepted sample types: serum, plasma
- B. All serum and plasma samples have been inactivated using TritonX-100 according to the following protocol prior to shipping to Olink:
 - 1. Maximum volume/sample for Triton inactivation: 90 µL
 - 2. Final TritonX-100 concentration: 1%
 - 3. Duration: 2 hours

I hereby certify that: (i) the COVID-19 samples provided to Olink for testing are provided in compliance with the above set forth requirements; and (ii) I am duly authorized to provide this Certification on behalf of my Institution.

By: _____

Name: _____

Title: _____

Institution name: _____

Date: _____