



Uppsala 2023-06-08

Change history

Olink® Explore Assay Nomenclature

| Date of Change | OlinkID | UniProt ID | Assay (=Gene Name) | Explore 384 Panel | Change |
|----------------|----------|------------|--------------------|--------------------|---|
| 2021-02-27 | OID20084 | O95684 | CEP43 | Cardiometabolic | Gene Name FGFR1OP changed to CEP43 |
| | OID20085 | Q12912 | IRAG2 | Cardiometabolic | Gene Name LRMP changed to IRAG2 |
| | OID20260 | P48960 | ADGRE5 | Cardiometabolic | Gene Name CD97 changed to ADGRE5 |
| | OID20903 | P14868 | DARS1 | Neurology | Gene Name DARS changed to DARS1 |
| | OID21044 | Q6P4E1 | GOLM2 | Neurology | Gene Name CASC4 changed to GOLM2 |
| | OID21068 | O76076 | CCN5 | Neurology | Gene Name WISP2 changed to CCN5 |
| | OID21086 | P12081 | HARS1 | Neurology | Gene Name HARS changed to HARS1 |
| | OID21209 | Q96NB1 | CEP20 | Oncology | Gene Name FOPNL changed to CEP20 |
| | OID21444 | O95388 | CCN4 | Oncology | Gene Name WISP1 changed to CCN4 |
| 2022-03-11 | OID30253 | PODPI2 | GATD3 | Cardiometabolic II | Gene Name GATD3B changed to GATD3 Uniprot ID A0A0B4J2D5 changed to PODPI2 |

www.olink.com

Olink Proteomics, Dag Hammarskjölds väg 52B

Uppsala Science Park, SE-751 83 Uppsala, Sweden

Phone: +46 (0)18 444 39 70, info@olink.com, Reg no: 559046-8632



| | | | | | |
|------------|----------|--------|-----------|-----------------|---|
| 2023-06-08 | OID20857 | Q9Y5P4 | CERT1 | Neurology | Gene Name CERT changed to CERT1 |
| | OID21084 | P23381 | WARS1 | Neurology | Gene Name WARS changed to WARS1 |
| | OID20125 | NA | NT-proBNP | Cardiometabolic | Gene Name NTproBNP change to NT-proBNP |